

SULFUR ANALYZER by U.V. FLUORESCENCE

ASTM D5453 / ASTM D6667 / ASTM D 7183 / ISO 20 846 NF M0759 - Licence TOTAL



PRINCIPLE

Sulfur Analysis

Mineralisation in gaseous phase of sulfur compounds forming SO₂ molecules detected by ultraviolet fluorescence (photomultiplier tube measures ultraviolet radiation emitted when excited SO₂ molecules return to base).

TECHNICAL CHARACTERISTICS

➤ Detector

- Graphic Liquid Crystal Display (LCD) / - Real-time synoptic flow diagram display
- Programmable measurement ranges / - Linearization of the measure
- Interactive menu driven software with enhanced speed display in 4 languages

➤ Apparatus

MINERALIZATION PART:

- A quartz combustion tube with Teflon output connector
- A combustion furnace at a constant controlled temperature of 1000°C (F2) for oxidation S → SO₂

A NO MEASUREMENT PART:

A specific detector of SO₂ by U.V. fluorescence.

SIGNAL/ CALCULATION/ STORAGE PROCESSING PART:

An associated computer which manages:

- SO₂ peaks integration / - Calculation of calibration coefficients / - Display of analysis results
- Storage on hard disk / - Automation and alarms.

GAS CONTROL PART:

Two gas circuits for Inert Gas and Oxygen, with gas pressure and flow regulators, pressure controllers and flowmeters.

ACCESSOIRES PART:

- A syringe-pusher for automatic injection of liquid specimens at controlled speed
- A colour printer for analysis and calibration results.

OPTIONS :

- Automatic injector for liquid samples
- Semiautomatic system for sampling and injection of gas in liquid LPG and / or gaseous phase

MINIMUM MAINTENANCE and easy to operate

VIDEO MONITOR: displays SO₂ peak as soon as the test starts

RESULTS (calibration or analysis) are automatically calculated at the end of the test and then printed

CALIBRATION "Monopoint" or Multipoint "(linear regression) mode selection by icon. The detector is linear; it is possible to dispense with the Multipoint calibration..

SWITCH FROM SOLID TO LIQUID MODE: in a few minutes.

APPLICATIONS

Liquid, solid or gaseous products, essentially petroleum products, but also organic compounds compatible with the method (industrial chemical products, rubber, synthetics, and so on).

TECHNICAL SPECIFICATIONS

DETECTION METHOD

Sulfer U.V. Fluorescence

SAMPLE SIZE

Liquid 20 to 100 µl with Syringe

Gas 1 to 25 ml with Syringe
10ml with Gaz / GPL Sampler

LPG 10 µl with sampler Gas /LPG

Solid 0.5 to 100 mg with boat

TYPICAL ANALYSIS TIME

Liquid and Gas About 5 minutes

Solid 5 to 10 minutes (depending upon type and weight of sample)

SUPPLY

Inert (Argon or Helium) 99.995% - 3 bar / 100 to 200 ml/min

Oxygen 99.998% - 3 bar / 200 to 300 ml/min

Electric 230 V – 50Hz – 1200 W

ACCURACY

At 0.5 ppm level +/- 0,05 ppm

At 1000 ppm level +/- 15 ppm

MEASUREMENT RANGE

DETECTION LIMITS


LOW QUANTIFICATION LIMIT

Sulfer	10 ppb to within approximately 10%	20 ppb
--------	---------------------------------------	--------

DIMENSIONS

Analyzer, Without computer 80 cm x 55 cm x 30 cm (BxHxD) / weight: 35 Kg

Schematic diagram of Gas circuits



We made also :



Wickbold / Sulfur

ASTM D2784 / ASTM D2785 -
AFNOR M41-009 / ISO 4260 -
IP243 / DIN 51408 NF. EN
24260 / EN41



***Tubular furnace with
temperature controller for
laboratory***



Chlorine Analyzer

AOX - Pox - Eox according to
ISO 9562



Nitrogen Analyzer

ASTM D4629 / ASTM D6069 /
ASTM D5176 / ASTM D7184 /
NF EN 12260 / NF M 07-058



Tri-four pyrolysis for

Tritium, Carbon 14, Chlorine 36,
Iodine 129

ERALY & Associés
4 Rue Georges Besse – Bât I
78330 FONTENAY-LE-FLEURY – France
Tel : +33(0)1 77 04 80 97
Fax : +33(0)1 77 04 80 96
www.eraly.com